 0373	Conforme alla direttiva Europea 93/42/CEE Complies with the European Directive 93/42 / EEC Cumple con la Directiva Europea 93/42/EEC Está en conformidad con la Directiva Europea 93/42/CEE Avrupa Direktifi 93/42/EEC ile uyumludur Comblinnann sé Treoir Eorpach 93/42/CEE Entspricht der europäischen Richtlinie 93/42/EWG		Non utilizzare se la confezione è danneggiata Do not use if package is damaged No utilizar si el embalaje está dañado. Nào utilize se a embalagem estiver danificada Amblajlaş baxarılmırsa durumda kullannmayın Nä húsidi ná dhéantar damáiste don phacáistiú Nicht verwenden, wenn die Verpackung beschädigt ist
	Sterile se non aperto o danneggiato / Raggi Beta Sterile if not opened or damaged / Beta rays Esteril si no está aberto o dañado / Beta Rays Esteril se não foi aberto ou estiver danificado / Raios Beta Aclannannas veyra hasa gáirmíse sterilid / Beta Imlarra Sterile mura n-oscailtear nó má dhéantar damáiste dó / Béite-ghathanna Steril, wenn ungeöffnet oder beschädigt / Betastrahlen		Riparare dalla pioggia e umidità Protect from rain and humidity Proteger de la lluvia y la humedad. Proteja da chuva e da humidade Yağmurdan ve nemden koruyun Cosain ó bhaisteach agus ó thaise Vor Regen und Feuchtigkeit schützen
	Non riutilizzare Do not reuse No reutilizar Nào reuse Tékar kullannayn Nä athúsaid Nicht wiederverwenden		Non resterilizzare Do not resterilize No reesterilizar Nào reesterilize Tékar sterilizet emeyin Nä athsterilíú Nicht erneut sterilisieren
	Consultare le istruzioni per l'uso Consult the instructions for use Consulte las instrucciones de uso Consulte as instruções de uso Kullann taliriallanna bakan Téigh i gcomhairle leis na treochra úsáide Beachten Sie die Gebrauchsanweisung		Contiene materiale biologico di origine animale Contains biological material of animal origin Contiene material biológico de origen animal. Contém material biológico de origem animal Hayvan kökenli biyolojik materyal içerir Tá ábhar bithéolaíoch de bhundú ainmleach ann Enthält biologisches Material tierischen Ursprungs
	Temperatura di conservazione Storage temperature Temperatura de almacenamiento Temperatura de armazenamento Depzallam siccadigh Teocht stóirála Lagertemperatur		Doppio confezionamento sterile Double sterile packaging Doble embalaje estéril Embalagem dupla estéril Cift steril amblaj Pacáistiú sterilid dubhalte Doppelt steril verpackt
	Numero di catalogo Catalog number Número de catálogo Número de catálogo Katalog númer Umíhr chatalóige Katalognummer		Tenere al riparo dalla luce solare Keep away from sunlight Mantener alejado de la luz solar Manter longe da luz solar Güneş ışığından uzak tutun Coinnigh ar sháid ó sholas na gréine Von Sonnenlicht fernhalten
	Codice partita / Numero lotto Match code / Lot number Código de lote / Número de lote Código do Lote / Número do Lote Parti kodu / Partı numarası Cód baice / Umíhr bhaice Chargencode / Chargennummer		Vettore identificazione dispositivo Device Identification Vector Vector de identificación del dispositivo Vetor de identificação do dispositivo Cihaz tanımlama vektörü Veicteoir aitheantais glais Geräteidentifikationsvektor
	Data di scadenza Expiration date Fecha de caducidad Data de validade Son kullanna tarihi Dáta éaga Verfallsdatum		
	Data di fabbricazione Manufacturing date Fecha de fabricación Data de fabricação Üretim tarihi Dáta déantaisíochta Herstellungsdatum		
	Dispositivo medico Medical device Dispositivo medico Dispositivo médico Tibbi cihaz Feiste leighis Medizinisches Gerät		

Fabbricante :
MAGGI srl Fraz. San pietro, 10
Piovà Massaia – 14026 (AT) Italy

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EXAFLEX

EXAFLEX High density collagen matrix, consisting of non crosslinked pericardium. Derived from bovines of maximum 20 months. The product allows the formation of newly formed collagen and the subsequent vascularization of the implant. The implant will subsequently be remodeled and replaced with host tissue. The product DOES NOT induce inflammation

INDICATIONS FOR USE:

The EXAFLEX membrane is intended for use as a prosthesis for the treatment of surgical repair of soft tissue deficiencies, positioning of breast implants. Breast reconstructive surgery.

CONTRAINDICATIONS:

Known or suspected sensitivity to derivatives of bovine origin. The product has not been tested on pregnant subjects. The product must not come into contact with the central nervous system or the eye. Do not use the product after the expiry date.

ADVERSE REACTIONS:

As with any surgical procedure, general risks include: infections, allergic reactions, edema, echimosis, foreign body reactions, acute or chronic inflammatory reactions, adhesions, seromas and relaxation of the repaired tissue. The patient must be aware of these risks and others associated with surgery in general and anesthesia.

NOTE TO THE DOCTOR:

The doctor must inform the patient of the indications, contraindications, warnings and precautions and possible complications indicated in this document.

PRECAUTIONS:

The product is disposable, therefore after opening the package any remaining products cannot be reused or resterilized. Before use, check the integrity of the packaging, otherwise do not use. Improper use could cause complications of a surgical nature. Do not use the product if not stored according to the manufacturer's instructions.

WARNING:

The product must never be dehydrated during use, use constant immersions in sterile physiological solution. In the event that the product must be used in contaminated sites, it is recommended to take appropriate local and systemic measures to contain contamination/infection.

STORAGE:

The product must be stored in clean and dry environments. The storage temperature must not be lower than 0° C / 32° F. DO NOT store in refrigerated cells.

TRANSPORT:

Transport in the intact package can take place under normal conditions. Both the product and the original packaging can be subjected to temperature changes for short periods.

INSTRUCTIONS:

The instructions below are not references for surgical technique and do not replace institutional protocols or professional clinical judgments regarding patient care. Verify that the product packaging is in good condition and that the packages containing the membrane are intact. Open the packages respecting the asepsis of the operating room. Rehydrate the membrane by immersion in sterile physiological solution at room temperature. The product can be soaked even for long periods, 2 or more hours. The product can be immersed in antibiotic-treated saline. Use the product according to surgical needs.

USE IN THE WAYS AND USES PERMITTED BY LOCAL REGULATIONS

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ÚSÁID AR NA HÚSÁIDÍ AGUS NA SEALAÍ A CEADAÍTEAR AG RIALACHÁIN ÁITIÚLA